







When all you need is everything; it is good to know there is one company that you can count on to deliver complete solutions and services. From Dental, Industry, Chemistry we develop better processes everyday - all backed by unrivaled customer service and support - we through our experience help you to raise game and increase your competitive edge.



the better way

All Resins Dimensionally Tested By Cube Printings from Corner to Center Of Printing Area

First test is dimensional stability test. We put 9 cubes as shown at the left picture to the printing table. Most of the printers has a light intensity difference from center to corner by %9-%30. We measure all the printed results before & after post cure and see the standart deviation. Thanks to our technology we get 9-16 micron standart deviation from the printed results. One of the most famous resin's standart deviation is more than 26 micron.

All resins needs to work as fast as possible

We control our resins formulations' reaction speed when the resin exposure by UV light. We achieve this not adding more photo initiators to the resin. By the way there is no color change after days and weeks. For example; with a light intensity 7 mW/cm^2; we are able work with 0,8-1,1 sec exposure times. A model which has around 300 layers (50 micron) can be printed in 19 minutes.

All resins needs to work by short separation distances not to loose much time

We develop our resins to work fast with all types of printers. To achieve that, we decrease the separation force of our resin from printing film. For example, one of the most famous printer brand works with its own resin by a 10 mm separation distance. We are able to work with 3.7 mm. It means %40-%50 less printing times.

All Resins Geometrically Tested By 3d Deformation Analysis

As a last testing; we make 3D deformation analysis with the original datas. Our model resins has below than 100 micron deformation for a full arch. This is a high level precision for the 3D printing market.

We Create our own formulations as to the process need

First we analyze the process and then start to develop formulations. We make our own synthesis for some of the indicators in our synthesis laboratory. Surface modifications are the key point to control the chemical bonds. We are professional regarding these critical points.

We make our in-house analysis and test as to the ISO standards

We measure all of our resins regarding the Elasticity modulus, flexural strength, tensile strength, hardness and elongation. Also water absorption tests, colour stability tests, phase separation tests are made in house from our analysis laboratory.

ArmaResin









UNDER THE BACK LIGHT

Opalessance Effect. Esthetic like ceramic materials.



This is just temp resin. Wait for the ArmaResin Permanent. Will be certified Q3 - 2023...

It will be again game changer...



6

"Biocompatible resin for dental models. Highly accurate and durable."



> Eliminated Phase Separation During Storage / Long Shelf Life

	Mechanical Properties			
Property	Procedure	Result		
Flexural Strength	ISO 4049	118,88		
Elasticity Modulus	ISO 4049	2047,05		
Viscosity (m.Pa.s 23 0C)		470		
Elongation at Break (mm)	ISO 4049	3,12		



- Standart Colour: Peach

- Alternative Colour: Gray



Optimised polymerisation during UV penetration Controlled Reaction Speed





Model HP is more than a model resin. It has scratch resistance and carvable like gypsum

Uniquely for all die models and implant models; this resin does not show any deformation. You are sure about the precision



For long time; this resin is not effected by sunlight. Even after months



Mechanical Properties

8

"Biocompatible resin for dental models. Highly accurate and cost effective."



> Eliminated Phase Separation During Storage / Long Shelf Life

Property	Procedure	Result			
Flexural Strength	ISO 4049	100,31			
Elasticity Modulus	ISO 4049	2081,3			
Viscosity (m.Pa.s 23 0C)		110			
Elongation at Break (mm)	ISO 4049	1,95			

Colour Alternatives

- Standart Colour: Beige

- Alternative Colour: Gray



Optimised polymerisation during UV penetration Controlled Reaction Speed





Economic resin for tooth restorations. Not only cheapest in the dental market; but also precise as well as the others.

It is dimensionally stable as well as the Model HP.



It has fast printing speeds. Dur to low viscosity and low separation forces; very easy to print with Model ST.



10

"Biocompatible resin for temporaies. Ultra high flexural strength. Low water absorption and stable colors."



> Eliminated Phase Separation During Storage / Long Shelf Life

	Me	chanical Prop
Property	Procedure	Resul
Flexural Strength	ISO 10477	125,46
Elasticity Modulus	ISO 10477	2329,48
Viscosity (m.Pa.s 23 0C)		400
Elongation at Break (mm)	ISO 10477	3,7
Water Absorption (mg/mm3)	ISO 4049	0,078
Colour Alternatives - A1-A2-A3-B1-BL-Bleach - Other colors are on demand		24 MONTHS
CF Medical Product Class I	On Progress	MA BPA

Controlled Reaction Speed







Reasonable price but outstanding properties. More than %95 correction for full arch without any support

Stabil colors before & after post cure. Ultra White - BL - B1 - A1 - A2 - A3 - A3,5 Opalessance Effect. Esthetic like glass ceramics



More than a temp resin. 125 Mpa flexural strength. Plus 4,5 mm elongation at break (ISO 10477)

A1 Colour After Post Cure

11



"Biocompatible resin for dentures. Ultra high flexural strength, low water absorption and stable colors."



Eliminated Phase Separation During Storage / Long Shelf Life

1 0 1			lechanical Prop
	Property	Procedure	Result
	Flexural Strength	ISO 10477	125,46
	Elasticity Modulus	ISO 10477	2329,48
	Viscosity (m.Pa.s 23 0C)		400
	Elongation at Break (mm)	ISO 10477	3,7
aResin	Water Absorption (mg/mm3)	ISO 4049	0,078
Polymer Resin 1.0 kg	Colour Alternatives - Original - Light - Other colors are on demand		24 MONTHS
S nm	CCE Medical Product Class Ila On Progress	On Progress	MMA BPA FREE FREE

Optimised polymerisation during UV penetration Controlled Reaction Speed

1000 gr



13



"Biocompatible resin for dental casting. Cost effective and fast working."



> Eliminated Phase Separation During Storage / Long Shelf Life

	Mechanical Properties			
Property	Procedure	Result		
Flexural Strength	ISO 4049	33,89		
Elasticity Modulus	ISO 4049	640,31		
Viscosity (m.Pa.s 23 0C)		300		
Elongation at Break (mm)	ISO 4049	7,11		
Burnout		< % 0,05		

Colour Alternatives



- Standart Colour: Transperent Yellowish

Optimised polymerisation during UV penetration Controlled Reaction Speed

1000 gr

Arma Resin

3D Printing Polymer Resin

DENTAL

CAST V 1.0

DLP/LCD 385-405 nm

1.0 ka

ARMA



Flexible before post-cure but dimensionally stable. Usable before & after post-cure

Able to work with the cheapest investments



"Biocompatible resin for surgical guides. Cost effective and highly precise"



Eliminated Phase Separation During Storage / Long Shelf Life

ISO 4049	88,17
ISO 4049	1937,54
·	300
ISO 4049	1,71
	24
On Progress	

Optimised polymerisation during UV penetration **Controlled Reaction Speed**



1000 gr

Transparent look Rigid material and robust during surgery

> Dimensionally stable No taste and neutral smell Autoclavable



"Biocompatible resin for occlusal splints. High flexural strength, No taste - Neutral Odour"



> Eliminated Phase Separation During Storage / Long Shelf Life

Property	Procedure	Result
Flexural Strength	ISO 4049	90
Elasticity Modulus	ISO 4049	2329,48
Viscosity (m.Pa.s 23 0C)		400
Elongation at Break (mm)	ISO 4049	4,5
Colour Alternatives		24
- Original - Transperent - Other colors are on demand		MONTHS

Optimised polymerisation during UV penetration Controlled Reaction Speed

1000 gr

Arma Resin

3D Printing Polymer Resin

SPLINT DLP/LCD 385-405 nm 1.0kg

CE





"Biocompatible resin for customized trays. Strong and durable to withstand the forces of taking the patient impression."

Specialities

- Ultra fast printing
- Eliminating stone models
- > High Flexural Strength.
- > Improves impression accuracy
- > Excellent Dimensional Stability

Mechanical Properties

- > Compatible with DLP / LCD
- > Works with 385 / 405 nm

> Eliminated Phase Separation During Storage / Long Shelf Life

Property	Procedure	Result
Flexural Strength	ISO 4049	143,03
Elasticity Modulus	ISO 4049	2640,76
Viscosity (m.Pa.s 23 0C)		470
Elongation at Break (mm)	ISO 4049	3,41



Colour Alternatives

- Standart Colour: Blue



Optimised polymerisation during UV penetration Controlled Reaction Speed



COST EFFECTIVE HIGH SPEED ACCURATE

COMPATIBLE WITH

SIGA®

SHINING 3D CREALITY

phrozen ANYCUBIC



MICRAFT



FLASHFORGE

and the other open system printers

ArmaResin











Arma Mach	🍪 Arma	Mold	H Arma	Resin	Arma	Shoe	ArmaAcryl
🔆 Arma	Chem	🏶 Arma	Jewel	🗰 ArmaC	Optics	Arma	Comp



GOSB Kemal Nehrozoğlu Cad. GOSB Teknoparkı 1. Üretim Binası 3-4-5. Birim Gebze - Kocaeli - Türkiye 41400

www.armadental.com

info@**armadental**.com

sinan.gunuc@armadental.com

Phone: +90 262 678 89 85 WP: +90 530 464 53 54